


Instruction sheet for Standard Microfluidic Mixing chip




Use:

This device serves as an on chip filtration device which may be used in conjunction with other Microfluidic devices

Materials provided:

Two pre assembled Standard Microfluidic Auto Dilution Chips including fittings.

Product Description:



The design for the mixer provides a junction for two solutions to combine along a path length of 12mm. The mixed solutions are filtered further down the micro channel through a restrictor and may be collected.

Note:

To place an order:

- By Phone: 1.877.776.8321
- By Fax: 304.293.4904
- Online: www.proteabio.com
orders@proteabio.com
- By Mail: Protea Biosciences, Inc.
955 Hartman Run Road #210
Morgantown, WV 26507